

aceRNA Technologies Co., Ltd.

CM detector RNA Switch™

PROTOCOL

Product name	Cat. No.	Size***	Storage
CM detector RNA Switch™	P-0002	9 µg	-80°C

*** 0.5 µg/µL

Additional materials required

- 24-well tissue culture plates
- 1.5 mL microcentrifuge tubes (RNase/DNase free, Sterile)
- Lipofectamine MessengerMAX™ Transfection Reagent (ThermoFisher)
- Reference iRFP670 mRNA - various fluorescents can be used except EGFP
mCherry mRNA (TriLink) is recommended
- Control detector RNA Switch™ (P-0001)

aceRNA Technologies Co., Ltd. (<http://acernatec.com>)

46-29 Yoshidashimoadachi-cho, Sakyo-ku, Kyoto, 606-8304, JAPAN

Technical support: info@acernatec.com

Cardiomyocyte detection using CM detector RNA Switch™

The following procedure is for the transfection of cell populations containing human cardiomyocytes with CM detector RNA Switch™ in one well of a 24-well plate. Transfection efficiency and mRNA expression can be influenced by the cell type being transfected.

RNA Switch™ is synthetic mRNA. Therefore, transfection should be performed in a RNase-free working-environment and mRNA should be diluted and aliquoted in RNase-free Water. Store them at -80°C.

1. Plate cell population containing human cardiomyocytes in a 24-well tissue culture plate at a density of $2.0 \sim 3.0 \times 10^5$ cells/well and incubate the cells at 37°C and 5% CO₂ 24 h prior to transfection.
2. 24 h after seeding, aspirate the medium and add 0.5 mL of fresh medium to wells and incubate the cells at 37°C and 5% CO₂ 1 to 2 hrs before transfection.
3. Thaw Control detector RNA Switch™ (P-0001), CM detector RNA Switch™, and reference iRFP670 mRNA on ice or at 4°C. Transfect the cells with 0.1 µg of Control detector RNA Switch™ or CM detector RNA Switch™ mixed with 0.1 µg of reference iRFP670 mRNA using 1 ~ 1.5 µL of Lipofectamine MessengerMAX™ Transfection Reagent (ThermoFisher) according to the manual. Incubate the cells at 37°C and 5% CO₂.
4. 4 h after transfection, discard the medium and add 0.5 mL of fresh medium. Incubate the cells at 37°C and 5% CO₂.
5. Next day, fluorescence expression can be imaged by image analyzer and harvested for FACS analysis.
 - Control detector RNA Switch™: EGFP
 - CM detector RNA Switch™: EGFP
 - Reference mRNA: iRFP670